

### 1-Step method (Reverse transfection method)

Timeline		Steps	Procedure Details									
Day 0	1	<p>Reagent Dilution Buffer</p> <p><b>Dilute ScreenFect™ mRNA Reagent*<sup>1</sup> in Dilution Buffer, and then mix well</b> *<sup>1</sup> Vortex the reagent before use</p>	Component		96-well		24-well		12-well		6-well	
	2	<p>mRNA Dilution Buffer</p> <p><b>Dilute mRNA in Dilution Buffer, and then mix well</b></p>	Dilution Buffer for ScreenFect™ mRNA		5 μL		25 μL		50 μL		125 μL	
	3	<p>mRNA-lipid complex</p> <p><b>Add diluted mRNA to diluted ScreenFect™ mRNA Reagent, and then incubate for 5 minutes ~ at room temperature*<sup>2</sup></b> *<sup>2</sup> Incubation is available until the step 4 has been completed</p>	mRNA : Transfection Reagent ratio		1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4
	4	<p>Cultured cells</p> <p><b>Prepare required cells for transfection</b></p>	ScreenFect™ mRNA Transfection Reagent		0.15 μL	0.2 μL	0.75 μL	1.0 μL	1.5 μL	2.0 μL	3.75 μL	5.0 μL
	5	<p>Cell suspension</p> <p><b>Detach cells and prepare the cell suspension, and then transfer the required numbers of cell suspension to cell culture plate</b></p>	Dilution Buffer for ScreenFect™ mRNA		5 μL		25 μL		50 μL		125 μL	
Day 1 ~	6	<p>mRNA-lipid complex</p> <p><b>Add mRNA-lipid complex from step 2 to well of cell culture plate from step 4</b></p>	mRNA (0.1-2.5 μg / μL)		50 ng		250 ng		500 ng		1250 ng	
		<p><b>Final composition [per well]</b></p>	Diluted mRNA		5 μL		25 μL		50 μL		125 μL	
		<p>Cell Detachment (Trypsin or Accutase®)</p>	Diluted ScreenFect™ mRNA Transfection Reagent		5 μL		25 μL		50 μL		125 μL	
		<p><b>Visualize/analyze transfected cells</b></p>	Adherent cells or suspension cells		1.0-4.0 × 10 <sup>4</sup>		0.5-2.0 × 10 <sup>5</sup>		1.0-4.0 × 10 <sup>5</sup>		0.25-1.0 × 10 <sup>6</sup>	
			mRNA amount		50 ng		250 ng		500 ng		1250 ng	
			ScreenFect™ mRNA Transfection Reagent used		0.15 or 0.2 μL		0.75 or 1.0 μL		1.5 or 2.0 μL		3.75 or 5.0 μL	
			Medium volume		100 μL		500 μL		1000 μL		2000 μL	
			Incubate cells for 1 day ~ at 37°C. Then, analyze transfected cells.									
			For support, please visit the <a href="http://screenfect.jp">http://screenfect.jp</a>									

### 2-Step method (Forward transfection method)

Timeline		Steps
1 Day 0		Seed cells to be 70-90% confluent at transfection
		<p>Dilute ScreenFect™ mRNA Reagent<sup>※1</sup> in Dilution Buffer, and then mix well  <sup>※1</sup> Vortex the reagent before use</p> <p>Dilute mRNA in Dilution Buffer, and then mix well</p>
3 Day 1		Add diluted mRNA to diluted ScreenFect™ mRNA Reagent, and then incubate for 5 minutes ~ at room temperature <sup>※2</sup> <sup>※2</sup> Incubation is available until the step 4 has been completed
		Add mRNA-lipid complex from step 3 to well of cell culture plate from step 1
5 Day 2 ~		Visualize/analyze transfected cells

Procedure Details		96-well	24-well	12-well	6-well				
Component									
Adherent cells or suspension cells		1.0-4.0 × 10 <sup>4</sup>	0.5-2.0 × 10 <sup>5</sup>	1.0-4.0 × 10 <sup>5</sup>	0.25-1.0 × 10 <sup>6</sup>				
Seed cells to be 70-90% confluent at transfection. The medium replacement may improve a transfection efficiency in a 2-Step method.									
Dilution Buffer for ScreenFect™ mRNA		5 μL	25 μL	50 μL	125 μL				
mRNA : Transfection Reagent ratio		1 : 3	1 : 4	1 : 3	1 : 4				
ScreenFect™ mRNA Transfection Reagent		0.15 μL	0.2 μL	0.75 μL	1.0 μL	1.5 μL	2.0 μL	3.75 μL	5.0 μL
Dilution Buffer for ScreenFect™ mRNA		5 μL	25 μL	50 μL	125 μL				
mRNA (0.1-2.5 μg / μL)		50 ng	250 ng	500 ng	1250 ng				
Diluted mRNA		5 μL	25 μL	50 μL	125 μL				
Diluted ScreenFect™ mRNA Transfection Reagent		5 μL	25 μL	50 μL	125 μL				
Final composition [per well]		96-well	24-well	12-well	6-well				
mRNA-lipid complex		10 μL	50 μL	100 μL	250 μL				
mRNA amount		50 ng	250 ng	500 ng	1250 ng				
ScreenFect™ mRNA Transfection Reagent used		0.15 or 0.2 μL	0.75 or 1.0 μL	1.5 or 2.0 μL	3.75 or 5.0 μL				
Medium volume		100 μL	500 μL	1000 μL	2000 μL				

Incubate cells for 1 day ~ at 37°C. Then, analyze transfected cells.