






### 1-Step method (Reverse transfection method)

Timeline	Steps
<b>1</b> 	<p><b>Dilute ScreenFect™ mRNA Reagent</b>※1 in Dilution Buffer, and then mix well            ※1 Vortex the reagent before use</p> <p><b>Dilute mRNA and SFA P-reagent</b> in Dilution Buffer, and then mix well</p>
	<p><b>2</b>   mRNA-lipid complex</p> <p><b>3</b>   Cultured cells</p> <p><b>4</b>   Cell suspension</p> <p><b>5</b>   mRNA-lipid complex</p>
<b>6</b> <b>Day 1 ~</b> 	<p><b>Visualize/analyze transfected cells</b></p>

Procedure Details		96-well		24-well		12-well		6-well	
Component									
Dilution Buffer for ScreenFect™ mRNA		5 μL		25 μL		50 μL		125 μL	
mRNA : Transfection Reagent ratio		1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4
ScreenFect™ mRNA Transfection Reagent		0.15 μL	0.2 μL	0.75 μL	1.0 μL	1.5 μL	2.0 μL	3.75 μL	5.0 μL
Dilution Buffer for ScreenFect™ mRNA		5 μL		25 μL		50 μL		125 μL	
mRNA (0.1-2.5 μg / μL)		50 ng		250 ng		500 ng		1250 ng	
SFA P-reagent (2μL / μg mRNA)		0.1 μL		0.5 μL		1.0 μL		2.5 μL	
Diluted mRNA (+SFA P-reagent)		5 μL		25 μL		50 μL		125 μL	
Diluted ScreenFect™ mRNA Transfection Reagent		5 μL		25 μL		50 μL		125 μL	
Adherent cells or suspension cells		1.0-4.0 × 10 <sup>4</sup>		0.5-2.0 × 10 <sup>5</sup>		1.0-4.0 × 10 <sup>5</sup>		0.25-1.0 × 10 <sup>6</sup>	
Cell Detachment (Trypsin or Accutase®)									
Final composition [per well]		96-well		24-well		12-well		6-well	
mRNA-lipid complex		10 μL		50 μL		100 μL		250 μL	
mRNA amount		50 ng		250 ng		500 ng		1250 ng	
SFA P-reagent used		0.1 μL		0.5 μL		1.0 μL		2.5 μL	
ScreenFect™ mRNA Transfection Reagent used		0.15 or 0.2 μL		0.75 or 1.0 μL		1.5 or 2.0 μL		3.75 or 5.0 μL	
Medium volume		100 μL		500 μL		1000 μL		2000 μL	
Incubate cells for 1 day ~ at 37°C. Then, analyze transfected cells.									

### 2-Step method (Forward transfection method)

Timeline		Steps	Procedure Details																																																														
1	Day 0	 Pre-Cultured cells	Seed cells to be 70-90% confluent at transfection	<table border="1"> <thead> <tr> <th>Component</th> <th colspan="2">96-well</th> <th colspan="2">24-well</th> <th colspan="2">12-well</th> <th colspan="2">6-well</th> </tr> </thead> <tbody> <tr> <td>Adherent cells or suspension cells</td> <td colspan="2">1.0-4.0 × 10<sup>4</sup></td> <td colspan="2">0.5-2.0 × 10<sup>5</sup></td> <td colspan="2">1.0-4.0 × 10<sup>5</sup></td> <td colspan="2">0.25-1.0 × 10<sup>6</sup></td> </tr> </tbody> </table>								Component	96-well		24-well		12-well		6-well		Adherent cells or suspension cells	1.0-4.0 × 10 <sup>4</sup>		0.5-2.0 × 10 <sup>5</sup>		1.0-4.0 × 10 <sup>5</sup>		0.25-1.0 × 10 <sup>6</sup>																																					
			Component	96-well		24-well		12-well		6-well																																																							
Adherent cells or suspension cells	1.0-4.0 × 10 <sup>4</sup>		0.5-2.0 × 10 <sup>5</sup>		1.0-4.0 × 10 <sup>5</sup>		0.25-1.0 × 10 <sup>6</sup>																																																										
Seed cells to be 70-90% confluent at transfection. The medium replacement may improve a transfection efficiency in a 2-Step method.																																																																	
2		 Reagent Dilution Buffer SFA P-reagent mRNA Dilution Buffer	Dilute ScreenFect™mRNA Reagent* <sup>1</sup> in Dilution Buffer, and then mix well ※1 Vortex the reagent before use  Dilute mRNA and SFA P-reagent in Dilution Buffer, and then mix well	<table border="1"> <tbody> <tr> <td>Dilution Buffer for ScreenFect™mRNA</td> <td colspan="2">5 μL</td> <td colspan="2">25 μL</td> <td colspan="2">50 μL</td> <td colspan="2">125 μL</td> </tr> <tr> <td>mRNA : Transfection Reagent ratio</td> <td>1 : 3</td> <td>1 : 4</td> <td>1 : 3</td> <td>1 : 4</td> <td>1 : 3</td> <td>1 : 4</td> <td>1 : 3</td> <td>1 : 4</td> </tr> <tr> <td>ScreenFect™mRNA Transfection Reagent</td> <td>0.15 μL</td> <td>0.2 μL</td> <td>0.75 μL</td> <td>1.0 μL</td> <td>1.5 μL</td> <td>2.0 μL</td> <td>3.75 μL</td> <td>5.0 μL</td> </tr> </tbody> </table>								Dilution Buffer for ScreenFect™mRNA	5 μL		25 μL		50 μL		125 μL		mRNA : Transfection Reagent ratio	1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4	ScreenFect™mRNA Transfection Reagent	0.15 μL	0.2 μL	0.75 μL	1.0 μL	1.5 μL	2.0 μL	3.75 μL	5.0 μL																											
		Dilution Buffer for ScreenFect™mRNA	5 μL		25 μL		50 μL		125 μL																																																								
mRNA : Transfection Reagent ratio	1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4	1 : 3	1 : 4																																																									
ScreenFect™mRNA Transfection Reagent	0.15 μL	0.2 μL	0.75 μL	1.0 μL	1.5 μL	2.0 μL	3.75 μL	5.0 μL																																																									
3	Day 1	 mRNA-lipid complex	Add diluted mRNA (+SFA P-reagent) to diluted ScreenFect™mRNA Reagent, and then incubate for 5 minutes ~ at room temperature* <sup>2</sup> ※2 Incubation is available until the step 4 has been completed	<table border="1"> <tbody> <tr> <td>Dilution Buffer for ScreenFect™mRNA</td> <td colspan="2">5 μL</td> <td colspan="2">25 μL</td> <td colspan="2">50 μL</td> <td colspan="2">125 μL</td> </tr> <tr> <td>mRNA (0.1-2.5 μg / μL)</td> <td colspan="2">50 ng</td> <td colspan="2">250 ng</td> <td colspan="2">500 ng</td> <td colspan="2">1250 ng</td> </tr> <tr> <td>SFA P-reagent (2μL / μg mRNA)</td> <td colspan="2">0.1 μL</td> <td colspan="2">0.5 μL</td> <td colspan="2">1.0 μL</td> <td colspan="2">2.5 μL</td> </tr> </tbody> </table>								Dilution Buffer for ScreenFect™mRNA	5 μL		25 μL		50 μL		125 μL		mRNA (0.1-2.5 μg / μL)	50 ng		250 ng		500 ng		1250 ng		SFA P-reagent (2μL / μg mRNA)	0.1 μL		0.5 μL		1.0 μL		2.5 μL																												
		Dilution Buffer for ScreenFect™mRNA	5 μL		25 μL		50 μL		125 μL																																																								
mRNA (0.1-2.5 μg / μL)	50 ng		250 ng		500 ng		1250 ng																																																										
SFA P-reagent (2μL / μg mRNA)	0.1 μL		0.5 μL		1.0 μL		2.5 μL																																																										
Add diluted mRNA (+SFA P-reagent) to diluted ScreenFect™mRNA Reagent, and then incubate for 5 minutes ~ at room temperature* <sup>2</sup> ※2 Incubation is available until the step 4 has been completed																																																																	
4		 mRNA-lipid complex	Add mRNA-lipid complex from step 3 to well of cell culture plate from step 1	<table border="1"> <tbody> <tr> <td>Diluted mRNA (+SFA P-reagent)</td> <td colspan="2">5 μL</td> <td colspan="2">25 μL</td> <td colspan="2">50 μL</td> <td colspan="2">125 μL</td> </tr> <tr> <td>Diluted ScreenFect™mRNA Transfection Reagent</td> <td colspan="2">5 μL</td> <td colspan="2">25 μL</td> <td colspan="2">50 μL</td> <td colspan="2">125 μL</td> </tr> </tbody> </table>								Diluted mRNA (+SFA P-reagent)	5 μL		25 μL		50 μL		125 μL		Diluted ScreenFect™mRNA Transfection Reagent	5 μL		25 μL		50 μL		125 μL																																					
		Diluted mRNA (+SFA P-reagent)	5 μL		25 μL		50 μL		125 μL																																																								
Diluted ScreenFect™mRNA Transfection Reagent	5 μL		25 μL		50 μL		125 μL																																																										
Add mRNA-lipid complex from step 3 to well of cell culture plate from step 1																																																																	
5	Day 2 ~		Visualize/analyze transfected cells	<table border="1"> <tbody> <tr> <td>Final composition [per well]</td> <td colspan="2">96-well</td> <td colspan="2">24-well</td> <td colspan="2">12-well</td> <td colspan="2">6-well</td> </tr> <tr> <td>mRNA-lipid complex</td> <td colspan="2">10 μL</td> <td colspan="2">50 μL</td> <td colspan="2">100 μL</td> <td colspan="2">250 μL</td> </tr> <tr> <td>mRNA amount</td> <td colspan="2">50 ng</td> <td colspan="2">250 ng</td> <td colspan="2">500 ng</td> <td colspan="2">1250 ng</td> </tr> <tr> <td>SFA P-reagent used</td> <td colspan="2">0.1 μL</td> <td colspan="2">0.5 μL</td> <td colspan="2">1.0 μL</td> <td colspan="2">2.5 μL</td> </tr> <tr> <td>ScreenFect™mRNA Transfection Reagent used</td> <td colspan="2">0.15 or 0.2 μL</td> <td colspan="2">0.75 or 1.0 μL</td> <td colspan="2">1.5 or 2.0 μL</td> <td colspan="2">3.75 or 5.0 μL</td> </tr> <tr> <td>Medium volume</td> <td colspan="2">100 μL</td> <td colspan="2">500 μL</td> <td colspan="2">1000 μL</td> <td colspan="2">2000 μL</td> </tr> </tbody> </table>								Final composition [per well]	96-well		24-well		12-well		6-well		mRNA-lipid complex	10 μL		50 μL		100 μL		250 μL		mRNA amount	50 ng		250 ng		500 ng		1250 ng		SFA P-reagent used	0.1 μL		0.5 μL		1.0 μL		2.5 μL		ScreenFect™mRNA Transfection Reagent used	0.15 or 0.2 μL		0.75 or 1.0 μL		1.5 or 2.0 μL		3.75 or 5.0 μL		Medium volume	100 μL		500 μL		1000 μL		2000 μL	
		Final composition [per well]	96-well		24-well		12-well		6-well																																																								
mRNA-lipid complex	10 μL		50 μL		100 μL		250 μL																																																										
mRNA amount	50 ng		250 ng		500 ng		1250 ng																																																										
SFA P-reagent used	0.1 μL		0.5 μL		1.0 μL		2.5 μL																																																										
ScreenFect™mRNA Transfection Reagent used	0.15 or 0.2 μL		0.75 or 1.0 μL		1.5 or 2.0 μL		3.75 or 5.0 μL																																																										
Medium volume	100 μL		500 μL		1000 μL		2000 μL																																																										
Incubate cells for 1 day ~ at 37°C. Then, analyze transfected cells.																																																																	