








1-Step method (Reverse transfection method)






Timeline	Steps
1  Reagent Dilution Buffer	Dilute ScreenFect™ A plus Reagent* ¹ in Dilution Buffer, and then mix well ※1 Vortex the reagent before use
	 siRNA Dilution Buffer
2  siRNA-lipid complex	Add diluted siRNA to diluted ScreenFect™ A plus Reagent, and then incubate for 5 minutes ~ at room temperature* ² ※2 Incubation is available until the step 4 has been completed
3  Cultured cells	Prepare required cells for transfection
4  Cell suspension	Detach cells and prepare the cell suspension, and then transfer the required numbers of cell suspension to cell culture plate
5  siRNA-lipid complex	Add siRNA-lipid complex from step 2 to well of cell culture plate from step 4
6 Day 1 ~ 	Visualize/analyze transfected cells

Procedure Details				
Component	96-well	24-well	12-well	6-well
Dilution Buffer for ScreenFect™ A plus	5 µL	25 µL	50 µL	125 µL
ScreenFect™ A plus Transfection Reagent	0.1 – 0.3 µL	0.5 – 1.5 µL	1.0 - 3.0 µL	2.5 – 7.5 µL
Dilution Buffer for ScreenFect™ A plus	5 µL	25 µL	50 µL	125 µL
siRNA	1 pmol siRNA	5 pmol siRNA	10 pmol siRNA	25 pmol siRNA
Diluted siRNA	5 µL	25 µL	50 µL	125 µL
Diluted ScreenFect™ A plus Transfection Reagent	5 µL	25 µL	50 µL	125 µL
Adherent cells or suspension cells	1.0-4.0 × 10 ⁴	0.5-2.0 × 10 ⁵	1.0-4.0 × 10 ⁵	0.25-1.0 × 10 ⁶
Cell Detachment (Trypsin or Accutase®)				
Final composition [per well]	96-well	24-well	12-well	6-well
siRNA-lipid complex	10 µL	50 µL	100 µL	250 µL
siRNA amount	1 pmol siRNA	5 pmol siRNA	10 pmol siRNA	25 pmol siRNA
ScreenFect™ A plus Transfection Reagent used	0.1 – 0.3 µL	0.5 – 1.5 µL	1.0 – 3.0 µL	2.5 – 7.5 µL
Medium volume	100 µL	500 µL	1000 µL	2000 µL
Incubate cells for 1 day ~ at 37°C. Then, analyze transfected cells.				
For support, please visit the http://screenfect.jp				

ScreenFect™ A plus Transfection Protocol

The optimal condition for successful transfection varies. Start any new transfection by testing the recommended two concentrations of ScreenFect™ A plus Reagent to determine an optimum amount.

2-Step method (Forward transfection method)

Timeline		Steps	Procedure Details					
1	Day 0	 <p>Pre-Cultured cells</p>	<p>Seed cells to be 70-90% confluent at transfection</p>	Component	96-well	24-well	12-well	6-well
				Adherent cells or suspension cells	1.0-4.0 × 10 ⁴	0.5-2.0 × 10 ⁵	1.0-4.0 × 10 ⁵	0.25-1.0 × 10 ⁶
Seed cells to be 70-90% confluent at transfection. The medium replacement may improve a transfection efficiency in a 2-Step method.								
2		 <p>Dilute ScreenFect™ A plus Reagent*¹ in Dilution Buffer, and then mix well ※¹ Vortex the reagent before use</p> <p>Dilute siRNA in Dilution Buffer, and then mix well</p>	Dilution Buffer for ScreenFect™ A plus	5 μL	25 μL	50 μL	125 μL	
			ScreenFect™ A plus Transfection Reagent	0.1 μL	0.5 μL	1.0 μL	2.5 μL	
3	Day 1	 <p>Add diluted siRNA to diluted ScreenFect™ A plus Reagent, and then incubate for 5 minutes ~ at room temperature*² ※² Incubation is available until the step 4 has been completed</p>	Dilution Buffer for ScreenFect™ A plus	5 μL	25 μL	50 μL	125 μL	
			siRNA	1 pmol siRNA	5 pmol siRNA	10 pmol siRNA	25 pmol siRNA	
4		 <p>Add siRNA-lipid complex from step 3 to well of cell culture plate from step 1</p>	Diluted siRNA	5 μL	25 μL	50 μL	125 μL	
			Diluted ScreenFect™ A plus Transfection Reagent	5 μL	25 μL	50 μL	125 μL	
Final composition [per well]				96-well	24-well	12-well	6-well	
siRNA-lipid complex				10 μL	50 μL	100 μL	250 μL	
siRNA amount				1 pmol siRNA	5 pmol siRNA	10 pmol siRNA	25 pmol siRNA	
ScreenFect™ A plus Transfection Reagent used				0.1 μL	0.5 μL	1.0 μL	2.5 μL	
Medium volume				100 μL	500 μL	1000 μL	2000 μL	
5	Day 2 ~	 <p>Visualize/analyze transfected cells</p>	Incubate cells for 1 day ~ at 37°C. Then, analyze transfected cells.					

For support, please visit the <http://screenfect.jp>